M^CKESSON

McKesson Wound*Care* Companion[™]

Product Resource Guide

uning uning

Information about our advanced wound care products and solutions.

McKesson Wound*Care* Companion[™]

This product resource guide has been developed to provide details about Advanced Wound Care products and information clinicians need to know in order to comply with CMS and OASIS requirements.

At McKesson Medical-Surgical, we understand that providing the best care is top of mind. That is why we have developed this comprehensive product guide to help you select the



right products for advanced wound and skin care. We know this can be an ongoing challenge and we are committed to finding ways to make the treatment of wounds less challenging for you and your clinical staff.

This resource guide offers you a broad choice between reputable, national brand suppliers as well as our own McKesson Brands product portfolio. It is conveniently divided by advanced skin and wound care subcategories and suppliers as well as provides product description, size and manufacturer item number.

In addition to this product resource guide, we also have other resources to assist clinicians with solutions for advanced skin and wound care.

Ask your McKesson Medical-Surgical Account Manager about taking advantage of these opportunities today or visit our full McKesson Wound*Care* Companion[™] program at **mms.mckesson.com/wound-care**

Guide Reference

ADVANCED SKIN CARE PRODUCT PORTFOLIO	
Advanced Skin Care	4
ADVANCED WOUND CARE PRODUCT PORTFOLIO	
Antifungals	8
Wound Cleansers and Measuring Guides	
Wound Cleansers	9
Wound Measuring Guides	10
No-Sting Liquid Film Skin Protectors	11
Adhesive Removers	12
Wound Dressings	
Collagens with or without Antimicrobial Silver	12
Xeroform	13
Hydrogels	13
Hydrophilic Wound Dressings	14
Transparent Films	14
Hydrocolloids	15
Adhesive and Non-Adhesive Foams with or without Antimicrobial Silver	16
Silicone Adhesive Foams with or without Antimicrobial Silver	17
Calcium Alginates and Hydrofibers without Antimicrobial Silver	18
Calcium Alginates and Hydrofibers with Antimicrobial Silver	19
MEDIHONEY [®]	20
Cover Dressings	
Bordered Gauze	21
Composite Dressings	22
Specialty Absorptive Dressings	22
Compression Dressings	
Unna Boots	23
Unna Boots with Calamine	23
Multi-Layer Compression Bandages	24
Wound Care Resources	25

Advanced Skin Care Product Portfolio

Maintaining skin integrity is the first step in preventing wounds. Properly selecting and using the right skin care product can be an ongoing challenge. That's why McKesson Medical-Surgical has made it easier for you to identify the proper products and is committed to finding ways to make the treatment of skin less challenging for you.

Our products are available in various units of measure - from single products to cases. Contact your McKesson Medical-Surgical Account Manager or check out mms.mckesson.com to order today.

ADVANCED SKIN CARE

THERA[™] is a specially formulated group of eight products that will help you address skin issues and protect^{*} and treat skin irritations.

INFUSED WITH SAFFLEXTM VITAMIN COMPOUND

Each product in the THERA[™] Advanced Skin Care line is enriched with an all-in-one proprietary blend of vitamins and nourishing ingredients called SAFFLEX[™]. It contains:

- Vitamin E acetate
- Vitamin C
- Vitamin B3
- Vitamin B5
- Vitamin B6
- Safflower seed oil
- Hyaluronic acid
- Lavender-Ylang fragrance
- Bisabolol (an extract of chamomile)
- Ginger root extract



*The skin protection claim does not apply to the Foaming Body Cleanser and the Antifungal Body Powder.



Advanced skin care

		E	BODY CLEANSER	S	BODY CREAM	BODY SHIELDS		BODY POWDER	
SE	CODUCT LECTION IART	FOAMING BODY CLEANSER	ANTIMICROBIAL BODY CLEANSER 0.13% Benzalkonium Chloride	MOISTURIZING BODY CLEANSER 1.75% Dimethicone	MOISTURIZING BODY CREAM 1.75% Dimethicone	MOISTURIZING BODY SHIELD 1.5% Dimethicone	DIMETHICONE BODY SHIELD 5.0% Dimethicone	CALAZINC BODY SHIELD 20.0% Calamine)	ANTIFUNGAL BODY POWDER 2.0% Miconazole Nitrate
	CLEANSES								
	No-rinse	•	•	•					
(0	Head-to-toe / Hair & body	•							
CLEANSES	Antiseptic / Reduces &		•						
AN	inhibits bacterial growth		•						
C.E.	Exfoliates	•							
	Peri wash	•	•	•		ļ			
	Moisturizes		•	•		ļ			
	Protects / Acts as a barrier			•					
MOISTURIZES	MOISTURIZES								
URI	Face			•	•	•			
DIST	Body			•	•	•	•		
ž	Long-lasting hydration			•	•	•	•		
	PROTECTS								
G	Dry skin			•	•	•	•		
Ë	Dry, cracked skin			•	•	•	•		
PROTECTS	Wet, cracked skin							•	
PRO	Itching				•	•	•	•	
	Against excess moisture			•		•	•	•	•
	Minor irritation				•	•	•	•	
E.	HELPS PREVENT								
HELPS PREVENT	Excess Trans Epidermal Water Loss (e-TEWL) by sealing in moisture			•	•	•	•		
ILPS P	Diaper rash caused by wetness from urine			•		•	•	•	
Ï	Diaper rash caused by wetness from stool					•	•	•	
	TREATS								
	Athlete's foot								•
TS	Jock itch		1					1	•
TREATS	Ringworm		1					İ	•
Ë	Fungal skin irritation		1					İ	•
	Diaper rash caused by		1			•	•	•	
	wetness and/or stool					-		-	
NO	SIZE	5 FL. OZ. (148 ML)	4 FL. OZ. (118 ML)	8 FL. OZ. (237 ML)	4 FL. OZ. (118 ML)	4 G PACKET	4 FL. OZ. (118 ML)	4 G PACKET	3 OZ. (85 G)
IAT	QTY	12/cs	24/cs	12/cs	12/cs	144/bx, 6 bx/cs	12/cs	144/bx, 6 bx/cs	12/cs
DRN	ITEM #	871345	871341	871343	871352	871349	871347	871351	871356
NFC	MFR #	116-BCLF5OZ	116-BCLA4OZ	116-BCLM8OZ	116-BCRM4OZ	116-BSM4G	116-BSD4OZ	116-BSC4G	116-BPA3OZ
ORDERING INFORMATION	SIZE	9 FL. OZ. (266 ML)	8 FL. OZ. (237 ML)		32 FL.OZ. (946 ML) WITH PUMP	4 OZ. (113 G)		4 OZ. (113 G)	
DEI	QTY	12/cs	12/cs		12/cs	12/cs		12/cs	
OR	ITEM #	871346	871342		871355	871348		871350	
	MFR #	116-BCLF9OZ	116-BCLA8OZ		116-BCRM32OZ	116-BSM4OZ		116-BSC4OZ	

ADVANCED SKIN CARE

ConvaTec offers a full line of soothing, high-quality skin care products, designed for routine preventive care as well as for compromised skin. Aloe Vesta[®] and Sensi-Care[®] products are designed to work together in a simple 1-2-3 system: Cleanse, Moisturize and Protect.

Aloe Vesta[®] products are designed for everyday, routine preventive skin care. They contain natural botanicals and are formulated to gently cleanse, moisturize and protect fragile skin. The products work together to help maintain skin integrity. They are excellent for fragile skin at the risk of breakdown as well as gentle and easy to use. Aloe Vesta[®] provides an effective barrier that seals out moisture, contains protein-rich emollients to moisturize and is non-sensitizing as well as fragrance free.

Sensi-Care[®] products help to manage skin compromised by exposure to urine and caustic stool residue. Effective but gentle for cleansing and protecting patients with sensitive and compromised skin. The products are formulated with less fragrance and added humectants to meet the needs of fragile skin, as well as provide extra protection for high-risk skin or weepy and denuded skin.



Advanced skin care

CLEANSE SIZE





4 oz.	Aloe Vesta [®] Body Wash & Shampoo	324604
8 oz.	Aloe Vesta [®] Body Wash & Shampoo	324609
1 L	Aloe Vesta [®] Body Wash & Shampoo, Bottle	324612
1 L	Aloe Vesta® Body Wash & Shampoo, Dispenser	122101
4 L	Aloe Vesta [®] Body Wash & Shampoo	324611
4 oz.	Sensi-Care® Body Wash & Shampoo	325304
8 oz.	Sensi-Care® Body Wash & Shampoo	325309
4 oz.	Aloe Vesta [®] Cleansing Foam	325204
8 oz.	Aloe Vesta [®] Cleansing Foam	325208
8 per pack	Aloe Vesta® Bathing Cloths	325521
4 oz.	Sensi-Care [®] Perineal/Skin Cleanser	324504
8 oz.	Sensi-Care® Perineal/Skin Cleanser	324509
4 oz.	Aloe Vesta® Perineal/Skin Cleanser	324704
8 oz.	Aloe Vesta [®] Perineal/Skin Cleanser	324709

MFR #

PRODUCT DESCRIPTION



MC	SIC	τu	RI	ZE

SIZE	PRODUCT DESCRIPTION	MFR #
3 oz.	Sensi-Care [®] Body Cream	324403
2 oz.	Aloe Vesta® Daily Moisturizer	324802
4 oz.	Aloe Vesta® Daily Moisturizer	324804
8 oz.	Aloe Vesta® Daily Moisturizer	324809

.....



PROTECT				
SIZE	PRODUCT DESCRIPTION	MFR #		
4 oz.	Sensi-Care® Protective Barrier	325614		
2 oz.	Aloe Vesta® Protective Ointment	324913		
8 oz.	Aloe Vesta® Protective Ointment	324908		
2.1 oz.	Aloe Vesta [®] Clear Barrier Spray	413401		

Advanced Wound Care Product Portfolio

Properly selecting and using the right wound care product can be an ongoing challenge. That's why McKesson Medical-Surgical has made it easier for you to identify the proper products and is committed to finding ways to make the treatment of wounds less challenging for you.

Our products are available in various units of measure - from single products to cases. Contact your McKesson Medical-Surgical Account Manager or check out mms.mckesson.com to order today.

ANTIFUNGALS Antifungals

Antifungals are agents that may be topically applied to the skin. They are available in powder, cream, lotion or ointment form that may also be formulated with ingredients that provide a moisture barrier. Indicated when a fungal or yeast infection is diagnosed or strongly suspected.

100	THERA™		
	SIZE	PRODUCT DESCRIPTION	MFR #
	3 oz.	THERA™ Antifungal Body Powder	116-BPA3OZ
	CONVATEC		

SIZE	PRODUCT DESCRIPTION	MFR #
2 oz.	Aloe Vesta [®] Clear Antifungal Ointment, Tube	325102
5 oz.	Aloe Vesta® Clear Antifungal Ointment, Tube	325105



SIZE	PRODUCT DESCRIPTION	MFR #
3 ¼ oz. Tube	Antifungal SECURA™ Cream (Extra thick Formula, Contains Zinc Oxide)	59432900

WOUND CLEANSERS AND MEASURING GUIDES

Wound Cleansers

Wound cleansers such as sterile saline or preserved solutions with surfactants or antimicrobials are used to clean wounds and remove debris and drainage from the wound surface.

	2
E.	
-	
Dermal Work Cleanser	-
19	

MCKESSON		
SIZE	PRODUCT DESCRIPTION	MFR #
4 oz.	Dermal Wound Cleanser	1718
8 oz.	Dermal Wound Cleanser	1719
16 oz.	Dermal Wound Cleanser	1720
7.1 oz.	Saline Wound Flush	37-6507



Anti-microbial

SIZE	PRODUCT DESCRIPTION	MFR #
4 oz.	Puracyn Plus™, Twist	6504
4 oz.	Puracyn Plus™, Pump	6505
8 oz.	Puracyn Plus™, Twist	6508
8 oz.	Puracyn Plus™, Pump	6509
16 oz.	Puracyn Plus™, Flip Spout	6516
16 oz.	Puracyn Plus™ OTC, Liquid Trigger	6517



CONVATEC

SIZE	PRODUCT DESCRIPTION	MFR #
20 ML	Shur-Clens [®] Skin Wound Cleansing Solution	121222
12 oz.	SAF-Clens [®] AF Dermal Wound Cleanser	159712



SIZE	PRODUCT DESCRIPTION	MFR #
8 oz.	Dermal Wound Cleanser	59449200
16 oz.	Dermal Wound Cleanser	449000

WOUND CLEANSERS AND MEASURING GUIDES

Wound Measuring Guides

MCKESSON

SIZE	IZE PRODUCT DESCRIPTION MFR		
Circular grids	 Wound Measuring Guides Disposable clear plastic, 5" x 7" sheet Circular bull's-eye grid measures wound diameter in centimeters Facilitates obtaining wound tracing Contains ample space to document name and date Designed for easy, single-use convenience 	533-30012100	
6" ruler	 Wound Measuring Guides Disposable white paper 6" ruler Calibrated in inches and centimeter equivalents Contains ample space to document name and date Designed for easy, single-use convenience 	533-LR/50	



	1 2 3 4 5 6 7 8	101 101 11	12 13 14 15 CENTIMETERS
	MSKESSON	PATIENT'S NAME	
Wound Measuring Guide		DATE	Rothined 3) Withease Medical Gargesal Inc. Rothinand, VA 20228 UR-4716-5213

NO-STING LIQUID FILM SKIN PROTECTORS NO-Sting Liquid Film Skin Protectors

Liquid skin protectors are formulations designed to protect vulnerable areas from the effects of 1) mechanical or chemical injury and 2) excessive moisture due to incontinence, perspiration or wound drainage. They form a transparent coating on the skin.



CO	NVA	TEC
	INVA	

CONTAILO		
SIZE	PRODUCT DESCRIPTION	MFR #
3 ML	Sensi-Care® Sting free Skin Barrier, Wipe	413501
1 ML	Sensi-Care® Sting free Skin Barrier, Foam Applicator	420793
1 ML	Sensi-Care [®] Sting free Skin Barrier, Foam Applicator	420794
3 ML	Sensi-Care® Sting free Barrier, Foam Applicator	420795
3 ML	Sensi-Care® Sting free Barrier, Foam Applicator	420796
28 ML	Sensi-Care® Sting free Skin Barrier, Spray, Bottle	420797
50 ML	Sensi-Care® Sting free Skin Barrier, Spray, Can	413502



SIZE	PRODUCT DESCRIPTION	MFR #
1 oz.	NO-STING SKIN-PREP™, Protective Dressing Spray, Liquid Film Forming Skin Protectant	66800709
3.0 ML	NO-STING SKIN-PREP™, Protective Dressing Swab, Liquid Film Forming Skin Protectant	59420700
1.0 ML	NO-STING SKIN-PREP™, Protective Dressing Wipe, Liquid Film Forming Skin Protectant	59420600

ADHESIVE REMOVER Adhesive Remover

Adhesive removers are solvents formulated to dissolve adhesives and assist in removing dressings, skin protectors, tapes and any product that adheres to the skin. Their composition permits gentle and painless cleansing of all types of adhesive residue. Indicated for removing hydrocolloid dressings, acrylic and rubber based adhesives. They may also be used to defat and cleanse the skin before applying an adhesive product.



CONVATEC

30 per boxSensi-Care® Sting Free Adhesive Remover Wipe41350050 MLSensi-Care® Sting Free Adhesive Releaser Spray413499150 MLSensi-Care® Sting Free Adhesive Releaser Spray420798	SIZE	PRODUCT DESCRIPTION	MFR #
	30 per box	Sensi-Care [®] Sting Free Adhesive Remover Wipe	413500
150 MI Sensi-Care® Sting Free Adhesive Beleaser Sprav 420798	50 ML	Sensi-Care® Sting Free Adhesive Releaser Spray	413499
	150 ML	Sensi-Care® Sting Free Adhesive Releaser Spray	420798



SMITH & NEPHEW

SIZE	PRODUCT DESCRIPTION	MFR #
50 per box	UNI-SOLVE™ Adhesive Remover Wipes	402300

WOUND DRESSINGS Collagens with or without Antimicrobial Silver

Collagen wound dressings are gels, pads, particles, pastes, powders, sheets or solutions; derived from bovine, porcine or avian sources. They encourage deposition and organization of newly formed collagen fibers and granulation tissue in the wound bed, stimulates new tissue development and wound debridement. Primary dressing for partial/full thickness, tunneling wounds, wounds with minimal to heavy drainage, skin grafts, donor sites, granulating or necrotic wounds. Benefits include absorbent, non-adherent, may be used in combination with topical agents, easy to apply and remove, moldable, may be used in infected wounds. However, collagens require secondary dressing and cannot be used for third degree burns or eschar.



wound dressings **Xeroform**

Primary dressing that contains 3% Bismuth Tribromophenate to provide bacteriostatic properties and help reduce risk of infection and deodorize the wound. Helps maintain a moist wound healing environment to assist in autolytic debridement for wounds with light or no drainage.

-	MCKESSON		
And Ling	SIZE	PRODUCT DESCRIPTION	MFR #
	1" x 8"	Xeroform Petrolatum Dressings	2201
	2" x 2"	Xeroform Petrolatum Dressings	2202
	4" × 4"	Xeroform Petrolatum Dressings	2206
	5" x 9"	Xeroform Petrolatum Dressings	2207

wound dressings Hydrogels

Used as a primary dressing for wounds with light or no drainage to add or maintain moisture at the wound bed and assist in autolytic debridement. To protect peri-wound skin from maceration, apply a moisture barrier ointment or barrier film. Available in amorphous, impregnated gauze and sheet forms.

1	MCKESSON		
1	SIZE	PRODUCT DESCRIPTION	MFR #
	1 oz.	Hydrogel Amorphous Wound Dressings	1721
Management Manual Deman	3 oz.	Hydrogel Amorphous Wound Dressings	1722
Tin	2" x 2"	Hydrogel Impregnated Gauze Dressings	61-53022
	4" × 4"	Hydrogel Impregnated Gauze Dressings	61-53044
	2" x 2"	Hydrogel Sheet Dressings	61-82022
	4" × 4"	Hydrogel Sheet Dressings	61-82044



CONVATEC

SIZE	PRODUCT DESCRIPTION	MFR #
3 oz.	SAF-Gel [®] Dermal Wound Gel with Alginate	145730



SMITH & NEPHEW				
SIZE	PRODUCT DESCRIPTION	MFR #		
15 g	INTRASITE™ Gel	66027311		
3 oz.	SOLOSITE™ Wound Gel	449600		
2" x 2"	SOLOSITE™ Conformable Wound Gel Dressing	59482300		
4" x 4"	SOLOSITE™ Conformable Wound Gel Dressing	59482400		

To order, visit mms.mckesson.com or call your McKesson Medical-Surgical Account Manager.

wound dressings Hydrophilic Wound Dressings

Hydrophilic Wound Dressing with Zinc, Vitamins A and B6, Calcium and Magnesium. Formulated with a patented zinc-nutrient formulation, and along with its balanced pH technology, provides a moist, mildly acidic environment conducive to wound healing.

 MCKESSON		
SIZE	PRODUCT DESCRIPTION	MFR #
3 oz.	Hydrophilic Wound Dressings	61-SPD03

wound dressings **Transparent Films**

Transparent Films are used as a primary or secondary dressing that helps maintain a moist wound surface, to assist in autolytic debridement in wounds with necrotic tissue. Thin and water-resistant, yet moisture vapor permeable, which aids in moisture management. This allows the wound site to be visually monitored.

in the second	146			
City of	a and the second	17		
	3.	관	123-	
	12	122	100	1
	-	. 17	-	i.
	-		2	
	6		12	
	_		1	

SIZE PR	ODUCT DESCRIPTION	MFR #
2 3/8" x 2 3/4" Trar	nsparent Film Dressings, Frame Delivery	4985
4" x 4 ¾" Trar	nsparent Film Dressings, Frame Delivery	4986
6" x 8" Trar	nsparent Film Dressings, Frame Delivery	4987



SIZE	PRODUCT DESCRIPTION	MFR #
2 3/8" x 2 3/4"	OPSITE FLEXIGRID™ Transparent Adhesive Film Dressing	66024628
4" x 4 ¾"	OPSITE™ Transparent Adhesive Film Dressing	4963
6" x 8"	OPSITE FLEXIGRID™ Transparent Adhesive Film Dressing	66024631

wound dressings Hydrocolloids

Used as a primary or secondary dressing for wounds with light to moderate drainage. Helps maintain a moist wound environment to assist in autolytic debridement. When in contact with wound drainage, the dressing forms a fluid/gel environment over the wound bed. Thin hydrocolloid is semi-transparent, allowing visualization of the wound.



MCKESSON

SIZE	PRODUCT DESCRIPTION	MFR #
2" x 2"	Hydrocolloid Dressings, Thin	1882
$4" \times 4"$	Hydrocolloid Dressings, Thin	1883
6" x 6"	Hydrocolloid Dressings, Thin	1884
6" x 7"	Hydrocolloid Dressings, Sacral	1885
2" x 2"	Hydrocolloid Dressings, Film Backing	1886
$4" \times 4"$	Hydrocolloid Dressings, Film Backing	1887
6" x 7"	Hydrocolloid Dressings, Film Backing	1888
$4" \times 4"$	Hydrocolloid Dressings, Foam Backing	1889
6" x 6"	Hydrocolloid Dressings, Foam Backing	1890



CONVATEC

SIZE	PRODUCT DESCRIPTION	MFR #
3" x 3"	DuoDERM [®] Hydrocolloid Dressing, Extra Thin	187901
4" x 4"	DuoDERM [®] Hydrocolloid Dressing, Extra Thin	187955
6" x 6"	DuoDERM [®] Hydrocolloid Dressing, Extra Thin	187957
4" × 4"	DuoDERM [®] Signal [®] Hydrocolloid Dressing	403326
6" x 7"	DuoDERM [®] Signal [®] Hydrocolloid Dressing	410501
4" x 4"	DuoDERM [®] CGF [®] Hydrocolloid Dressing	187660
6" x 6"	DuoDERM [®] CGF [®] Hydrocolloid Dressing	187661



SIZE	PRODUCT DESCRIPTION	MFR #
2" x 2 ¾"	REPLICARE™ Thin Hydrocolloid Dressing	59484000
3 ½" x 5 ½"	REPLICARE™ Thin Hydrocolloid Dressing	59484100
6" x 8"	REPLICARE™ Thin Hydrocolloid Dressing	59484200
4" × 4"	REPLICARE™ Ultra Hydrocolloid Dressing	59484600
7" x 8"	REPLICARE™ Ultra Hydrocolloid Dressing	59484900

WOUND DRESSINGS

Adhesive and Non-Adhesive Foams with or without Antimicrobial Silver

Used as a primary or secondary dressing for wounds with moderate to heavy drainage. Highly absorbent central hydrocellular foam layer effectively absorbs wound drainage. Helps maintain a moist healing environment at the wound surface.



MCKESSON

SIZE	PRODUCT DESCRIPTION	MFR #
3" x 3"	Hydrocellular Foam Dressings, Acrylic Adhesive	4670
4" x 4"	Hydrocellular Foam Dressings, Acrylic Adhesive	4671
6" x 6"	Hydrocellular Foam Dressings, Acrylic Adhesive	4672
7" x 7"	Hydrocellular Foam Dressings, Acrylic Adhesive	4673
4" x 4"	Hydrocellular Foam Dressings, Non-Adhesive, Film Backing	4737
6" x 6"	Hydrocellular Foam Dressings, Non-Adhesive, Film Backing	4738
4" × 4"	Hydrophilic Foam Dressings, Non-Adhesive	4739



CONVATEC

SIZE	PRODUCT DESCRIPTION	MFR #
4" × 4"	AQUACEL [®] Foam, Non-Adhesive	420633
6" x 6"	AQUACEL [®] Foam, Non-Adhesive	420635
6" x 8"	AQUACEL [®] Foam, Non-Adhesive	420637
2" x 2"	AQUACEL [®] Ag Foam, Non-Adhesive	420639
4" x 4"	AQUACEL [®] Ag Foam, Non-Adhesive	420642
6" x 6"	AQUACEL [®] Ag Foam, Non-Adhesive	420645
8" x 8"	AQUACEL [®] Ag Foam, Non-Adhesive	420646
6" x 8"	AQUACEL® Ag Foam, Non-Adhesive	420806



U		
SIZE	PRODUCT DESCRIPTION	MFR #
3" x 3"	ALLEVYN™ Adhesive Hydrocellular Foam Dressing	66020043
7" x 7"	ALLEVYN™ Adhesive Hydrocellular Foam Dressing	66020045
4" x 4"	ALLEVYN™ Non-Adhesive Hydrocellular Foam Dressing	66927637
6" x 6"	ALLEVYN™ Non-Adhesive Hydrocellular Foam Dressing	66020093

WOUND DRESSINGS

Silicone Adhesive Foams with or without Antimicrobial Silver

Highly conformable primary or secondary dressing for wounds with moderate drainage. Perforated, soft silicone gel adhesive provides gentle, but firm, adhesion. Breathable, yet moisture and bacteria-resistant outer film layer seals the wound margins to help prevent drainage while absorbent inner hydrocellular foam pad maintains a moist wound environment to help promote better healing.



MCKESSON

SIZE	PRODUCT DESCRIPTION	MFR #			
3" x 3"	Hydrocellular Foam Dressings, Silicone Adhesive	4842			
$4" \times 4"$	Hydrocellular Foam Dressings, Silicone Adhesive	4843			
6" x 6"	Hydrocellular Foam Dressings, Silicone Adhesive	4844			
7" x 7"	Hydrocellular Foam Dressings, Silicone Adhesive, Sacral	4845			



CONVATEC SIZE PRODUCT DESCRIPTION MFR # 3" x 3" AQUACEL[®] Foam 420804 AQUACEL[®] Foam 4" x 4" 420680 5" x 5" AQUACEL[®] Foam 420619 7" x 7" AQUACEL® Foam 420621 8" x 5.5" AQUACEL® Foam, Heel 420625 8" x 7" AQUACEL[®] Foam. Small Sacral 420626 9.4" x 8.4" AQUACEL[®] Foam, Large Sacral 420828 3" x 3" AQUACEL[®] Ag Foam 420805 4" x 4" AQUACEL[®] Ag Foam 420681 5" x 5" AQUACEL® Ag Foam 420627 7" x 7" AQUACEL® Ag Foam 420628 8" x 5.5" AQUACEL® Ag Foam, Heel 420647 AQUACEL® Ag Foam, Small Sacral 8" x 7" 420648 9.4" x 8.4" AQUACEL® Ag Foam, Large Sacral 420827



SMITH & NEPHEW

SIZE	PRODUCT DESCRIPTION	MFR #
3" x 3"	ALLEVYN™ Gentle Border Hydrocellular Foam Dressing	66800276
4" × 4"	ALLEVYN™ Gentle Border Hydrocellular Foam Dressing	66800270
5" x 5"	ALLEVYN™ Gentle Border Hydrocellular Foam Dressing	66800279
6 5⁄/8" x 6 3⁄/4"	ALLEVYN™ Gentle Border Sacrum Hydrocellular Foam Dressing	66800898
5" X 5"	ALLEVYN™ Ag Gentle Adhesive Border	66800453
3" X 3"	ALLEVYN™ Ag Gentle Adhesive Border	66800452

To order, visit mms.mckesson.com or call your McKesson Medical-Surgical Account Manager.

wound dressings **Calcium Alginates and Hydrofibers** without Antimicrobial Silver

Primary dressing that can be cut to fit wounds with moderate to heavy drainage. They absorb wound fluid while maintaining a moist wound environment. The alginate, derived from seaweed, and the hydrofibers both turn to a moist gel mass material when saturated. Requires a secondary dressing.



MCKESSON		
SIZE	PRODUCT DESCRIPTION	MFR #
³ ⁄4" × 12"	Calcium Alginate Dressings, Rope	3564
2" x 2"	Calcium Alginate Dressings, Sheet	3561
4" × 4 ³ /4"	Calcium Alginate Dressings, Sheet	3562
4" x 8"	Calcium Alginate Dressings, Sheet	3563

AQUACEL® Hydrofiber® Wound Dressing

AQUACEL[®] EXTRA[™] Hydrofiber[®] dressing

AQUACEL[®] EXTRA™ Hydrofiber[®] dressing

MFR #

420127

420671

420674

PRODUCT DESCRIPTION





•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•

CONVATEC SIZE

³⁄₄" x 12"

2" x 2"

4 ¹/₄" x 4 ¹/₄"

SMITH & NEPHEW					
SIZE	PRODUCT DESCRIPTION	MFR #			
³ ⁄4" x 12"	ALGISITE™ M Calcium Alginate Dressing	59480400			
2" x 2"	ALGISITE™ M Calcium Alginate Dressing	59480100			
4" × 4"	ALGISITE™ M Calcium Alginate Dressing	59480200			
6" x 8"	ALGISITE™ M Calcium Alginate Dressing	59480300			

WOUND DRESSINGS Calcium Alginates and Hydrofibers with Antimicrobial Silver

Primary dressing that can be cut to fit wounds with moderate to heavy drainage. Assists with infection reduction and decreasing bioburden. Silver ions are released in the presence of wound fluid and protect the dressing from bacterial colonization and provide an effective barrier to bacterial penetration. Requires a secondary dressing.



MCKESSON

SIZE	PRODUCT DESCRIPTION	MFR #
³ ⁄4" x 12"	Calcium Alginate Dressings, Rope, with Antimicrobial Silver	3560
2" x 2"	Calcium Alginate Dressings, Sheet, with Antimicrobial Silver	3557
4" x 4 ¾"	Calcium Alginate Dressings, Sheet, with Antimicrobial Silver	3558
4" x 8"	Calcium Alginate Dressings, Sheet, with Antimicrobial Silver	3559



CONVATEC	
SIZE	PRODUCT DESCRIPTION
2" x 2"	AQUACEL® Ag EXTRA

2" x 2"	AQUACEL® Ag EXTRA	420675
4" × 5"	AQUACEL® Ag EXTRA	420677
6" x 6"	AQUACEL® Ag EXTRA	420678
8" x 12"	AQUACEL® Ag EXTRA	420679
³ ⁄4" x 18"	AQUACEL® Ag Ribbon	403771

MFR #



SMITH & NEPHEW

SIZE	PRODUCT DESCRIPTION	MFR #
³ ⁄4" x 12"	ACTICOAT™ Absorbent Silver-Coated Antimicrobial Calcium Alginate Dressing	20181
4" x 5"	ACTICOAT™ Absorbent Silver-Coated Antimicrobial Calcium Alginate Dressing	20381

WOUND DRESSINGS **MEDIHONEY**®

If debridement and the promotion of healthy tissue is your clinical objective, consider MEDIHONEY® for your wound. MEDIHONEY®'s high osmotic effect pulls fluid from deeper tissues to the wound surface promoting autolytic debridement. The naturally low pH of MEDIHONEY® helps lower the pH of the wound environment, which can impact chronic and stalled wounds. To protect peri-wound skin from maceration, apply a moisture barrier ointment or barrier film.

MEDIHONEY[®] **Gel:** Primary dressing for wounds with light to moderate drainage.

MEDIHONEY[®] **Honeycolloids**[™]**:** Primary dressing for wounds with light to moderate drainage.

MEDIHONEY® Calcium Alginates: Primary dressing for wounds with moderate to heavy drainage. Can be cut to fit within the wound.

MEDIHONEY®

SIZE	PRODUCT DESCRIPTION	MFR #
0.5 oz.	MEDIHONEY® Gel Tube Twist Off	31805
1.5 oz.	MEDIHONEY® Gel Tube Twist Off	31815
3 ½" x 3 ½"	MEDIHONEY® Honeycolloid™ Dressing, Adhesive	31422
6" × 6"	MEDIHONEY [®] Honeycolloid™ Dressing, Adhesive	31445
2" x 2"	MEDIHONEY [®] Honeycolloid™ Dressing, Non-adhesive	31222
4" x 5"	MEDIHONEY [®] Honeycolloid™ Dressing, Non-adhesive	31245
³ ⁄4" x 12"	MEDIHONEY® Calcium Alginate Rope Dressing	31012
2" x 2"	MEDIHONEY® Calcium Alginate Sheet Dressing	31022
4" x 5"	MEDIHONEY® Calcium Alginate Sheet Dressing	31045



MEDIHONEY[®] is a registered trademark of Medihoney Pty Ltd. and licensed to Derma Sciences, Inc. MEDIHONEY[®] is not a private label product of McKesson but is available for distribution through McKesson's exclusive brand of Advanced Wound Care products.

cover dressings Bordered Gauze

Primary or secondary bordered non-adherent gauze dressing for wounds with light to moderate drainage. One-step application is an excellent alternative to gauze and tape. Non-adherent central gauze pad provides for absorption of light to moderate drainage.



MCKESSON					
SIZE	PRODUCT DESCRIPTION	MFR #			
4" × 4"	Island Dressings, Non-Sterile	16-89244			
6" x 6"	Island Dressings, Non-Sterile	16-89266			
2" x 2"	Island Dressings, Sterile	16-89022			
4" × 4"	Island Dressings, Sterile	16-89044			
4" × 6"	Island Dressings, Sterile	16-89046			
4" × 8"	Island Dressings, Sterile	16-89048			
4" × 10"	Island Dressings, Sterile	16-89041			
4" × 14"	Island Dressings, Sterile	16-89042			
6" x 6"	Island Dressings, Sterile	16-89066			
6" x 8"	Island Dressings, Sterile	16-89068			



SIZE	PRODUCT DESCRIPTION	MFR #
4" × 4"	COVRSITE™ Cover Dressing	59714100
6" x 6"	COVRSITE™ Cover Dressing	59714400

cover dressings Composite Dressings

Primary or secondary bordered, non-adherent, occlusive gauze composite dressing for wounds with light to moderate drainage. Soft, flexible and comfortable for enhanced wearability. Water-resistant film layer acts as effective barrier to outside contamination and permits showering. Offers the convenience of an all-in-one dressing by combining an absorbent non-adherent central gauze pad surrounded by a breathable adhesive border, with a polyurethane film barrier.





MCKESSON

SIZE	PRODUCT DESCRIPTION	MFR #
4" x 4"	Barrier Island Dressings, Water-Resistant, Sterile	16-89144
6" x 6"	Barrier Island Dressings, Water-Resistant, Sterile	16-89166

SMITH & NEPHEW

SIZE	PRODUCT DESCRIPTION	MFR #
$4" \times 4"$	COVRSITE™ PLUS Composite Dressing	59715000
6" x 6"	COVRSITE™ PLUS Composite Dressing	59715100

cover dressings Specialty Absorptive Dressings

Can be used as a primary or secondary dressing for moderate to heavy drainage. Highly absorbent dressings help to decrease the risk of maceration and keep the harmful wound components away from the wound and peri-wound skin. Contains protective outer layer and non-adherent wound contact layer for easy removal.





MCKESSON

SIZE	PRODUCT DESCRIPTION	MFR #
4" x 5"	Super Absorbent Dressings, Sterile	61-89545
6" x 9"	Super Absorbent Dressings, Sterile	61-89569

SIZE	PRODUCT DESCRIPTION	MFR #
4" x 6	" EXU-DRY™ Dressing	5999004120
6" x 9	" EXU-DRY™ Dressing	5999006

COMPRESSION DRESSINGS Unna Boots

Gauze bandage that has been evenly impregnated with a non-hardening paste of zinc oxide to provide venous ulcer compression therapy. Both Unna Boots and Unna Boots with Calamine feature non-raveling, thread-locked gauze helping prevent loose threads that can contaminate wounds. Contains no gelatin.



Μ	CKESSON	

SIZE	PRODUCT DESCRIPTION	MFR #
3"	Unna Boot with Zinc Oxide, Non-Sterile	2066
4"	Unna Boot with Zinc Oxide, Non-Sterile	2067



CONVATEC

SIZE	PRODUCT DESCRIPTION	MFR #
3"	Unna-Flex [®] Unna Boot	650940
4"	Unna-Flex [®] Unna Boot	650941



SMITH & NEPHEW

SIZE	PRODUCT DESCRIPTION	MFR #		
3"	VISCOPASTE™ PB7 Zinc Paste Bandage	4956		

COMPRESSION DRESSINGS Unna Boots with Calamine

Gauze bandage that has been evenly impregnated with a non-hardening paste of zinc oxide and calamine to provide venous ulcer compression therapy and help ease itching.

	-
•	Unna Boot 7"
ł	100 000
Į	0.00
ſ	F139+12.3



MCKESSON

SIZE	PRODUCT DESCRIPTION	MFR #
3"	Unna Boot with Zinc Oxide & Calamine, Non- Sterile	2068
4"	Unna Boot with Zinc Oxide & Calamine, Non- Sterile	2069

COMPRESSION DRESSINGS Multi-Layer Compression Bandage Systems

Multi-Layer Compression Bandage Systems are superior bandaging systems for optimal compression and absorbency in managing venous ulcer therapy. Multi-layered system simultaneously absorbs and protects, while providing graduated sustained compression.



MCKESSON

MOREGOON			
SIZE	PRODUCT DESCRIPTION	MFR #	
N/A	Compression Bandage System, 2 Layer, Non-Sterile	1006	
N/A	Compression Bandage System, 4 Layer, Non-Sterile	1008	



SMITH & NEPHEW

SIZE	PRODUCT DESCRIPTION	MFR #
N/A	PROFORE™ Lite Multi-Layer Reduced Compression Bandaging System	66000771
N/A	PROFORE™ Multi-Layer Compression Bandaging System, 4 layer	66020016

Our products are available in various units of measure - from single products to cases. Contact your Account Manager for more information.

woundcare resources WoundCare Resources

In addition to the McKesson Wound*Care* Companion[™] program, please reference the following resources to assist in wound care education, prevention and treatment:



Association for the Advancement of Wound Care *aawconline.org*



McKesson Education Portal mms.mckesson.com/ online-education

Panel

www.npuap.org

Smith & Nephew

www.smith-nephew.com

National Pressure Ulcer Advisory



Center for Medicare & Medicaid Services www.CMS.gov



ConvaTec www.convatec.com







Wound Care Education Institute *www.wcei.net*



Kestrel Wound Source www.woundsource.com



McKesson Brands mms.mckesson.com/ mckesson-brand



> smith&nephew

Wound, Ostomy, Continence Nurses Society www.wocn.org

Notes

.....

Notes

.....

MSKESSON

Have clinical questions or need educational support?

McKesson Medical-Surgical's Board-Certified Wound, Ostomy, Continence Nurses are ready to help you today. **Call our clinicians at 877.611.0081.**

To learn more about WoundCare Companion™

visit mms.mckesson.com/wound-care

or call your McKesson Medical-Surgical Account Manager.

Be advised that information contained herein is intended to serve as a useful reference for informational purposes only and is not complete clinical information. This information is intended for use only by competent healthcare professionals exercising judgment in providing care. McKesson Medical-Surgical cannot be held responsible for the continued currency of or for any errors or omissions in the information. If any condition noted in these informational guidelines persists or worsens, reassess condition and determine whether further or different treatment is required.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers, not McKesson Medical-Surgical or its affiliates ("McKesson"), and have not been independently verified by McKesson. McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use. Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

McKesson Medical-Surgical is providing HCPCS codes based upon information provided by the product manufacturers. McKesson Medical-Surgical does not make representations about the accuracy of this information. The HCPCS codes do not indicate coverage or payment under any federal or state healthcare program. The entity that submits claims to any federal or state healthcare program has the responsibility to verify the accuracy of HCPCS codes it submits for payment purposes.

mms.mckesson.com/wound-care

All trademarks and registered trademarks are the property of their respective owners. © 2015 McKesson Medical-Surgical Inc. 2015-0125